

Certificate of Analysis

Product: Microalbumin Control
 Catalogue Number: BXC328B
 Batch Number: 162781
 Manufactured Date: 2016/11
 Expiry: 2017/12

This product has been tested by Quality Control and has met all final QC specifications. This product is released for packing and is suitable for use until the stated expiry date.

This product was produced according to ISO 13485 & ISO 9001 guidelines.

Lot	Method	Target	Range	Units
Level 1	Immunoturbidimetry	44.4	37.7 - 51.1	mg/L
Level 2	Immunoturbidimetry	88.7	75.3 - 102	mg/L

Products containing Human Source Material have been tested at the donor level for HbsAg and antibodies to HIV1, HIV2 and HCV and found to be non-reactive.

Signed: *Cherith Morris* (QC Manager)

Date: 2016/11