



**CERTIFICATE OF ANALYSIS**  
**FINISHED KIT QUALITY CONTROL**

PRODUCT NAME: **ALBUMIN** MANUFACTURE DATE: 03/10/2011  
 CATALOG #: A502-480  
 LOT #: 31052 EXP. DATE: 02/2013

COMPONENTS: ALBUMIN REAGENT LOT#: D110101011 EXP. DATE: 01/2014  
 ALBUMIN STANDARD LOT#: 24180 EXP. DATE: 02/2013

BLANK: 0.441 AT 600 NM

**TOLERANCE LEVELS OF SERUM CONTROLS**

<u>SERA ID</u>	<u>EXP. DATE</u>	<u>TOLERANCE RANGE</u> (On TC-84)	<u>TARGET</u> (g/dl)
15281	09/2012	3.7-4.6	4.1
15272	03/2012	1.9 - 2.9	2.4

**TEST VALUES/RESULTS**

<u>SERA ID</u>	<u>1</u>	<u>2</u>	<u>MEAN</u>	<u># PASS</u>
15281	4.0	3.9	4.0	2/2
15272	2.2	2.3	2.3	2/2

COMMENTS: VISUAL CHECK OK.

APPROVED BY:  DATE: 03/10/2011  
 BY: Quality Control