

TECO DIAGNOSTICS



1268 N. Lakeview Ave., Anaheim, CA 92807, U.S.A.
 Phone: (714) 463-1111 Toll Free (USA): (800) 222-9880
 Fax: (714) 463-1169
www.tecodiag.com

CERTIFICATE OF ANALYSIS
FINISHED KIT QUALITY CONTROL

PRODUCT NAME: ALBUMIN	MANUFACTURE DATE: 07/02/2015
CATALOG #: A502-480	
LOT #: 61765	EXP. DATE: 02/2017

COMPONENTS:	ALBUMIN REAGENT	LOT#: 52776	EXP. DATE: 02/2017
	ALBUMIN STANDARD	LOT#: 52777	EXP. DATE: 02/2017

TOLERANCE LEVELS OF SERUM CONTROLS

<u>SERA ID</u>	<u>EXP. DATE</u>	<u>TOLERANCE RANGE</u> (On TC-Matrix)	<u>TARGET</u> (g/dl)
27621	12/2015	3.9 - 4.8	4.3
27622	12/2015	2.5 - 3.1	2.8

TEST VALUES/RESULTS

<u>SERA ID</u>	<u>MEAN</u>	<u>PASS/FAIL</u>
27621	4.6	PASS
27622	3.0	PASS

COMMENTS: VISUAL CHECK OK.

APPROVED BY:

Jacqueline Sanchez
 Quality Control

DATE:

07/06/2015