



CERTIFICATE OF ANALYSIS
(FINISHED KIT QUALITY CONTROL)

PRODUCT NAME: **ALKALINE PHOSPHATASE (KINETIC)** MANUFACTURE DATE: 01/11/2010
 CATALOG #: A505-240G
 LOT #: 22257 EXP. DATE: 07/2011

COMPONENT: ALKALINE PHOSPHATASE REAGENT LOT#: 17706 EXP. DATE: 07/2011

REAGENT BLANK VS. WATER: 0.358 AT 405 NM

TOLERANCE LEVELS OF SERUM CONTROLS

<u>SERA ID</u>	<u>EXP. DATE</u>	<u>TOLERANCE RANGE</u> (On TC-84)	<u>TARGET</u> (IU/L)
15271	03/2012	59.3 – 110.2	85.0
15272	03/2012	252 – 360	468

TEST VALUES/RESULTS

<u>SERA ID</u>	<u>1</u>	<u>2</u>	<u>MEAN</u>	<u>#PASS</u>
15271	85	90.4	87.7	2/2
15272	396.7	335.8	366.3	2/2

COMMENTS: VISUAL CHECK OK.

APPROVED BY:	 Quality Control	DATE:	01/11/2010
-----------------	---------------------	-------	------------