



**CERTIFICATE OF ANALYSIS**  
**(FINISHED KIT QUALITY CONTROL)**

PRODUCT NAME: **ALKALINE PHOSPHATASE (KINETIC)** MANUFACTURE DATE: 02/05/2010  
 CATALOG #: A505-240  
 LOT #: 24130 EXP. DATE: 08/2012

COMPONENT: ALKALINE PHOSPHATASE REAGENT LOT#: 23319 EXP. DATE: 08/2012

REAGENT BLANK VS. WATER: 0.346 AT 405 NM

**TOLERANCE LEVELS OF SERUM CONTROLS**

<u>SERA ID</u>	<u>EXP. DATE</u>	<u>TOLERANCE RANGE</u> (On TC-84)	<u>TARGET</u> (IU/L)
15271	03/2012	59.3-110.2	85.0
15272	03/2012	252-468	360

**TEST VALUES/RESULTS**

<u>SERA ID</u>	<u>1</u>	<u>2</u>	<u>MEAN</u>	<u>#PASS</u>
15271	90.3	86.3	88.3	2/2
15272	387.3	363.4	375.4	2/2

COMMENTS: VISUAL CHECK OK.

APPROVED  
 BY:

Quality Control

DATE: 03/30/2010