



CERTIFICATE OF ANALYSIS
(FINISHED KIT QUALITY CONTROL)

PRODUCT NAME: **ALKALINE PHOSPHATASE (KINETIC)** MANUFACTURE DATE: 05/18/10
CATALOG #: A505-240
LOT #: 25580 EXP. DATE: 10/12

COMPONENT: ALKALINE PHOSPHATASE REAGENT LOT#: 25309 EXP. DATE: 10/12

REAGENT BLANK VS. WATER: AT 405 NM

TOLERANCE LEVELS OF SERUM CONTROLS

<u>SERA ID</u>	<u>EXP. DATE</u>	<u>TOLERANCE RANGE</u> (On TC-84)	<u>TARGET</u> (IU/L)
15271	03/2012	59.3-110.2	85.0
15272	03/2012	252-468	360

TEST VALUES/RESULTS

<u>SERA ID</u>	<u>1</u>	<u>2</u>	<u>MEAN</u>	<u>#PASS</u>
15271	85.1	72.1	78.6	2/2
15272	366.2	352.2	359.2	2/2

COMMENTS: VISUAL CHECK OK.

APPROVED
BY:


Quality Control

DATE: 05/18/10