



CERTIFICATE OF ANALYSIS
(FINISHED KIT QUALITY CONTROL)

PRODUCT NAME: ALKALINE PHOSPHATASE (KINETIC) MANUFACTURE DATE: 03/10/2011
CATALOG #: A505-240
LOT #: 30344 EXP. DATE: 04/2013

COMPONENT: ALKALINE PHOSPHATASE REAGENT LOT#: 28466 EXP. DATE: 04/2013

REAGENT BLANK VS. WATER: 0.016 AT 405 NM

TOLERANCE LEVELS OF SERUM CONTROLS

<u>SERA ID</u>	<u>EXP. DATE</u>	<u>TOLERANCE RANGE</u> (On TC-84)	<u>TARGET</u> (IU/L)
15281	09/2012	85.6-128	107
15272	03/2012	252-468	360

TEST VALUES/RESULTS

<u>SERA ID</u>	<u>1</u>	<u>2</u>	<u>MEAN</u>	<u>#PASS</u>
15281	122.4	119.8	121.1	2/2
15272	373.5	360.1	366.8	2/2

COMMENTS: VISUAL CHECK OK.

APPROVED
BY:

Quality Control

DATE: 03/11/2011