



CERTIFICATE OF ANALYSIS
(FINISHED KIT QUALITY CONTROL)

PRODUCT NAME: **ALT (SGPT) - KINETIC** MANUFACTURE DATE: 04/23/2010
CATALOG #: A525-240
LOT#: 24700 EXP. DATE: 04/2013

COMPONENT: ALT (SGPT) REAGENT LOT#: 24922 EXP. DATE: 04/2013

REAGENT BLANK VS. WATER: 1.850 AT 340 NM

TOLERANCE LEVELS OF SERUM CONTROLS

<u>SERA ID</u>	<u>EXP. DATE</u>	<u>TOLERANCE RANGE</u>	<u>TARGET (IU/L)</u>
15271	3/2012	24.8 – 38.0	31
15272	3/2012	68.4 – 102.6	85.5

TEST VALUES/RESULTS

<u>SERA ID</u>	<u>1</u>	<u>2</u>	<u>MEAN</u>	<u>#PASS</u>
15271	31.1	30.0	30.6	2/2
15272	91.9	90.9	90.4	2/2

COMMENTS: VISUAL CHECK OK.

APPROVED
BY:


Quality Control

DATE: 04/23/2010