



**CERTIFICATE OF ANALYSIS**  
**(FINISHED KIT QUALITY CONTROL)**

PRODUCT NAME: **ALT (SGPT) - KINETIC**      MANUFACTURE DATE: 10/25/2010  
CATALOG #: A525-400  
LOT#: 28477      EXP. DATE: 08/2013

COMPONENT: ALT (SGPT) REAGENT    LOT#: 26300    EXP. DATE: 08/2013

REAGENT BLANK VS. WATER: 1.819    AT 340 NM

**TOLERANCE LEVELS OF SERUM CONTROLS**

<u>SERA ID</u>	<u>EXP. DATE</u>	<u>TOLERANCE RANGE</u>	<u>TARGET</u> (IU/L)
15271	03/2012	24.8 – 38.0	31.0
15272	03/2012	69.3 – 103.8	86.6

**TEST VALUES/RESULTS**

<u>SERA ID</u>	<u>1</u>	<u>2</u>	<u>MEAN</u>	<u>#PASS</u>
15271	34.5	31.7	33.1	2/2
15272	87.4	91.1	89.2	2/2

COMMENTS: VISUAL CHECK OK.

APPROVED  
BY:

  
Quality Control

DATE: 10/25/2010