



CERTIFICATE OF ANALYSIS
(FINISHED KIT QUALITY CONTROL)

PRODUCT NAME: **ALT (SGPT) - KINETIC** MANUFACTURE DATE: 02/07/2011
 CATALOG #: A525-400
 LOT#: 29812 EXP. DATE: 11/2013

COMPONENT: ALT (SGPT) REAGENT LOT#: 29357 EXP. DATE: 11/2013

REAGENT BLANK VS. WATER: 1.784 AT 340 NM

TOLERANCE LEVELS OF SERUM CONTROLS

<u>SERA ID</u>	<u>EXP. DATE</u>	<u>TOLERANCE RANGE</u>	<u>TARGET</u> (IU/L)
15281	09/2012	28-42	35.0
15272	03/2012	69.3 – 103.8	86.6

TEST VALUES/RESULTS

<u>SERA ID</u>	<u>1</u>	<u>2</u>	<u>MEAN</u>	<u>#PASS</u>
15281	37.6	33.3	35.45	2/2
15272	85.6	87.3	86.45	2/2

COMMENTS: VISUAL CHECK OK.

APPROVED
BY:


Quality Control

DATE: 02/07/2011