

TECO DIAGNOSTICS



1268 N. Lakeview Ave., Anaheim, CA 92807, U.S.A.
Phone: (714) 463-1111 Toll Free (USA): (800) 222-9880
Fax: (714) 463-1169
www.tecodiag.com

CERTIFICATE OF ANALYSIS
(FINISHED KIT QUALITY CONTROL)

PRODUCT NAME: ALT (SGPT) - KINETIC MANUFACTURE DATE: 03/05/2014
CATALOG #: A525-400
LOT#: 53526 EXP. DATE: 12/2015

COMPONENT: ALT (SGPT) REAGENT LOT#: 43633 EXP. DATE: 12/2015

TOLERANCE LEVELS OF SERUM CONTROLS

| <u>SERA ID</u> | <u>EXP. DATE</u> | <u>TOLERANCE RANGE</u> (on TC Matrix) | <u>TARGET</u> (IU/L) |
|----------------|------------------|--|-------------------------|
| 27611 | 05/2015 | 26.0 – 39.0 | 32.5 |
| 27602 | 03/2015 | 71.6 – 107.4 | 89.5 |

TEST VALUES/RESULTS

| <u>SERA ID</u> | <u>MEAN</u> | <u>PASS/FAIL</u> |
|----------------|-------------|------------------|
| 27611 | 35 | PASS |
| 27602 | 93 | PASS |

COMMENTS: VISUAL CHECK OK.

APPROVED BY: Susan
Susan Tang
Quality Control

DATE: 03/06/2014