

TECO DIAGNOSTICS



1268 N. Lakeview Ave., Anaheim, CA 92807, U.S.A.
Phone: (714) 463-1111 Toll Free (USA): (800) 222-9880
Fax: (714) 463-1169
www.tecodiag.com

CERTIFICATE OF ANALYSIS
(FINISHED KIT QUALITY CONTROL)

PRODUCT NAME: ALT (SGPT) - KINETIC MANUFACTURE DATE: 07/06/2015
CATALOG #: A525-400
LOT#: 61717 EXP. DATE: 12/2017

COMPONENT: ALT (SGPT) REAGENT LOT#: 57705 EXP. DATE: 12/2017

TOLERANCE LEVELS OF SERUM CONTROLS

<u>SERA ID</u>	<u>EXP. DATE</u>	<u>TOLERANCE RANGE</u> (on TC 84)	<u>TARGET</u> (IU/L)
27621	12/2015	20.6 – 36.8	28.7
27622	12/2015	71.2 – 106.8	89.0

TEST VALUES/RESULTS

<u>SERA ID</u>	<u>MEAN</u>	<u>PASS/FAIL</u>
27621	32.9	PASS
27622	90.0	PASS

COMMENTS: VISUAL CHECK OK.

APPROVED BY:


Jacqueline Sanchez
Quality Control

DATE:

