



**CERTIFICATE OF ANALYSIS**  
**(FINISHED KIT QUALITY CONTROL)**

PRODUCT NAME: **AMYLASE - KINETIC**      MANUFACTURE DATE: 06/11/2010  
CATALOG #: A533-60  
LOT #: 26108      EXP. DATE: 03/2013

COMPONENT: AMYLASE REAGENT      LOT#: 24449      EXP. DATE: 03/2013

REAGENT BLANK VS. WATER: 0.083      AT 405 NM

**TOLERANCE LEVELS OF SERUM CONTROLS**

<u>SERA ID</u>	<u>EXP. DATE</u>	<u>TOLERANCE RANGE</u>	<u>TARGET</u> <u>(IU/L)</u>
15271	03/2012	51.7 – 63.2	57.0
15272	03/2012	268.5 – 328.1	298.0

**TEST VALUES/RESULTS**

<u>SERA ID</u>	<u>1</u>	<u>2</u>	<u>MEAN</u>	<u># PASS</u>
15271	58.9	58.3	58.6	2/2
15272	300.3	292.5	296.4	2/2

COMMENTS: VISUAL CHECK OK.

APPROVED  
BY:

  
Quality Control

DATE : 06/11/2010