



CERTIFICATE OF ANALYSIS
(FINISHED KIT QUALITY CONTROL)

PRODUCT NAME: **AMYLASE - KINETIC** MANUFACTURE DATE: 11/16/2010
CATALOG #: A533-60
LOT #: 28474 EXP. DATE: 11/2013

COMPONENT: AMYLASE REAGENT LOT#: 28552 EXP. DATE: 11/2013

REAGENT BLANK VS. WATER: 0.057 AT 405 NM

TOLERANCE LEVELS OF SERUM CONTROLS

| <u>SERA ID</u> | <u>EXP. DATE</u> | <u>TOLERANCE RANGE</u> | <u>TARGET</u> <u>(IU/L)</u> |
|----------------|------------------|------------------------|--------------------------------|
| 15271 | 03/2012 | 51.7 – 63.2 | 57.0 |
| 15272 | 03/2012 | 268.5 – 328.1 | 298.0 |

TEST VALUES/RESULTS

| <u>SERA ID</u> | <u>1</u> | <u>2</u> | <u>MEAN</u> | <u># PASS</u> |
|----------------|----------|----------|-------------|---------------|
| 15271 | 57.8 | 60.8 | 59.3 | 2/2 |
| 15272 | 321.5 | 322.1 | 321.8 | 2/2 |

COMMENTS: VISUAL CHECK OK.

APPROVED BY:  DATE: 11/16/2010
BY: Quality Control :