



**CERTIFICATE OF ANALYSIS**  
**(FINISHED KIT QUALITY CONTROL)**

PRODUCT NAME: **AMYLASE - KINETIC**      MANUFACTURE DATE: 12/06/2011  
 CATALOG #: A533-60  
 LOT #: 35945      EXP. DATE: 12/2014

COMPONENT: AMYLASE REAGENT      LOT#: 35921      EXP. DATE: 12/2014

REAGENT BLANK VS. WATER: 0.039      AT 405 NM

**TOLERANCE LEVELS OF SERUM CONTROLS**


<u>SERA ID</u>	<u>EXP. DATE</u>	<u>TOLERANCE RANGE</u>	<u>TARGET</u> (IU/L)
15281	09/2012	53.6 – 80.4	67.0
15282	09/2012	242.2 – 363.3	302.8

**TEST VALUES/RESULTS**

<u>SERA ID</u>	<u>1</u>	<u>2</u>	<u>MEAN</u>	<u># PASS</u>
15281	60.7	58.6	59.65	2/2
15282	303.9	307.2	305.55	2/2

COMMENTS: VISUAL CHECK OK.

APPROVED BY:

  
 Nour Qasqas  
 Quality Control Lab. Director

DATE: 12/06/2011