

TECO DIAGNOSTICS



1268 N. Lakeview Ave., Anaheim, CA 92807, U.S.A.
Phone: (714) 463-1111 Toll Free (USA): (800) 222-9880
Fax: (714) 463-1169
www.tecodiag.com

CERTIFICATE OF ANALYSIS
(FINISHED KIT QUALITY CONTROL)

PRODUCT NAME: **AST (SGOT) KINETIC** MANUFACTURE DATE: 07/08/2010
CATALOG #: A560-240
LOT #: 26582 EXP. DATE: 07/2012

COMPONENT: **AST (SGOT) REAGENT** LOT #: 19939 EXP. DATE: 07/2012

REAGENT BLANK VS. WATER: 1.656 AT 340 NM


TOLERANCE LEVELS OF SERUM CONTROLS

<u>SERA ID</u>	<u>EXP. DATE</u>	<u>TOLERANCE RANGE</u>	<u>TARGET (IU/L)</u>
15271	03/2012	26.4 – 52.0	39.0
15272	03/2012	152.4 – 228.6	190.5

TEST VALUES/RESULTS

<u>SERA ID</u>	<u>1</u>	<u>2</u>	<u>MEAN</u>	<u># PASS</u>
15271	44.7	32.1	38.4	2/2
15272	178.0	186.0	182.0	2/2

COMMENTS: VISUAL CHECK OK.

APPROVED BY:	 Quality Control	DATE:	<u>07/09/2010</u>
--------------	--	-------	-------------------