



CERTIFICATE OF ANALYSIS
(FINISHED KIT QUALITY CONTROL)

PRODUCT NAME: **AST (SGOT) KINETIC** MANUFACTURE DATE: 01/25/2011
 CATALOG #: A560-240
 LOT #: 28375 EXP. DATE: 11/2013

COMPONENT: AST (SGOT) REAGENT LOT #: 28256 EXP. DATE: 11/2013

REAGENT BLANK VS. WATER: 1.842 AT 340 NM

TOLERANCE LEVELS OF SERUM CONTROLS

<u>SERA ID</u>	<u>EXP. DATE</u>	<u>TOLERANCE RANGE</u>	<u>TARGET (IU/L)</u>
15271	03/2012	26.4 – 52.0	39.0
15272	03/2012	152.4 – 228.6	190.5

TEST VALUES/RESULTS

<u>SERA ID</u>	<u>1</u>	<u>2</u>	<u>MEAN</u>	<u># PASS</u>
15271	36.2	38.5	37.35	2/2
15272	166.0	173.1	169.55	2/2

COMMENTS: VISUAL CHECK OK.

APPROVED BY:	 Quality Control	DATE:	<u>01/25/2011</u>
--------------	--	-------	-------------------