



CERTIFICATE OF ANALYSIS
(FINISHED KIT QUALITY CONTROL)

PRODUCT NAME: **AST (SGOT) KINETIC** MANUFACTURE DATE: 11/30/2010
CATALOG #: A560-240
LOT #: 28887 EXP. DATE: 11/2013

COMPONENT: AST (SGOT) REAGENT LOT #: 28256 EXP. DATE: 11/2013

REAGENT BLANK VS. WATER: 1.740 AT 340 NM


TOLERANCE LEVELS OF SERUM CONTROLS

<u>SERA ID</u>	<u>EXP. DATE</u>	<u>TOLERANCE RANGE</u>	<u>TARGET (IU/L)</u>
15271	03/2012	26.4 – 52.0	39.0
15272	03/2012	152.4 – 228.6	190.5

TEST VALUES/RESULTS

<u>SERA ID</u>	<u>1</u>	<u>2</u>	<u>MEAN</u>	<u># PASS</u>
15271	32.8	31.8	32.3	2/2
15272	172.4	169.5	170.95	2/2

COMMENTS: VISUAL CHECK OK.

APPROVED BY:	 Quality Control	DATE:	<u>12/01/2010</u>
--------------	--	-------	-------------------