

TECO DIAGNOSTICS



1268 N. Lakeview Ave., Anaheim, CA 92807, U.S.A.
Phone: (714) 463-1111 Toll Free (USA): (800) 222-9880
Fax: (714) 463-1169
www.tecodiag.com

CERTIFICATE OF ANALYSIS
(FINISHED KIT QUALITY CONTROL)

PRODUCT NAME: AST (SGOT) MANUFACTURE DATE: 08/27/2014
CATALOG #: A560-400
LOT #: 56674 EXP. DATE: 11/2015

COMPONENT: AST (SGOT) REAGENT LOT #: 43636 EXP. DATE: 11/2015

TOLERANCE LEVELS OF SERUM CONTROLS


<u>SERA ID</u>	<u>EXP. DATE</u>	<u>TOLERANCE RANGE</u> (on TC-Matrix)	<u>TARGET</u> (U/L)
27611	05/2015	34.4 – 51.6	43.0
27602	03/2015	160.4 – 240.6	200.5

TEST VALUES/RESULTS

<u>SERA ID</u>	<u>MEAN</u>	<u>PASS/FAIL</u>
27611	42.0	PASS
27602	191.7	PASS

COMMENTS: VISUAL CHECK OK.

APPROVED BY:


Jacqueline Sanchez
Quality Control

DATE:

8/27/2014