



CERTIFICATE OF ANALYSIS
(FINISHED KIT QUALITY CONTROL)

PRODUCT NAME: **BILIRUBIN DIRECT** MANUFACTURE DATE: 08/01/2014
 CATALOG #: B538-480
 LOT #: 56462 EXP. DATE: 06/2017

COMPONENTS: 1. DIRECT BILIRUBIN REAGENT LOT#: 55353 EXP. DATE: 06/2017
 2. NITRITE REAGENT LOT#: 56009 EXP. DATE: 07/2017
 3. CALIBRATOR (5 mg/dL) LOT#: 56010 EXP. DATE: 07/2017

TOLERANCE LEVELS OF SERUM CONTROLS


<u>SERA ID</u>	<u>EXP DATE</u>	<u>TOLERANCE RANGE</u>	<u>TARGET</u> (mg/dl)
27611	05/2015	0.2-0.5	0.3
27602	03/2015	0.8-1.6	1.2

TEST VALUES/RESULTS

<u>SERA ID</u>	<u>MEAN</u>	<u>PASS /FAIL</u>
27611	0.3	PASS
27602	1.2	PASS

COMMENTS: VISUAL CHECK OK.

APPROVED BY:


 Jacqueline Sanchez
 Quality Control

DATE: 08/04/2014