



CERTIFICATE OF ANALYSIS
(FINISHED KIT QUALITY CONTROL)

PRODUCT NAME: **BILIRUBIN TOTAL** MANUFACTURE DATE: 06/27/2011
 CATALOG #: B576-480
 LOT #: 33198 EXP. DATE: 05/2014

<u>COMPONENTS:</u>	<u>CATALOG #</u>	<u>LOT#</u>	<u>EXP DATE</u>
1. TOTAL BILIRUBIN REAGENT	B576-A	32773	06/2014
2. NITRITE REAGENT	B538-B	33169	07/2014
3. CALIBRATOR (5.0 mg/dl)	B538-C	32109	05/2014

BLANK: -0.005 AT 550 nm

TOLERANCE LEVELS OF SERUM CONTROLS

<u>SERA ID</u>	<u>EXP DATE</u>	<u>TOLERANCE RANGE</u> (On TC-84)	<u>TARGET</u> (mg/dl)
15281	09/2012	1.0-1.7	1.2
15282	09/2012	4.9-7.3	6.1

TEST VALUES/RESULTS

<u>SERA ID</u>	<u>1</u>	<u>2</u>	<u>MEAN</u>	<u># PASS</u>
15281	1.1	1.1	1.1	2/2
15282	4.9	4.9	4.9	2/2

COMMENTS: VISUAL CHECK OK.

APPROVED BY:  DATE: 06/28/2011