



CERTIFICATE OF ANALYSIS
(FINISHED KIT QUALITY CONTROL)

PRODUCT NAME: **UREA NITROGEN (BUN) KINETIC** MANUFACTURE DATE: 01/08/2009
 CATALOG #: B550-400
 LOT #: 17427 EXP. DATE: 02/2011

COMPONENTS: 1. BUN REAGENT LOT#: 16302 EXP. DATE: 02/2011
 2. BUN STANDARD LOT#: 16710 EXP. DATE: 03/2011

REAGENT BLANK VS. WATER: 1.781 AT 340 NM

TOLERANCE LEVELS OF CONTROLS

<u>SERA ID</u>	<u>EXP. DATE</u>	<u>TOLERANCE RANGE</u> (On TC-84)	<u>TARGET</u> (mg/dl)
15221	08/2009	13.4 – 16.8	15.1
15242	08/2010	35.2 – 44.8	40.0

TEST VALUES/RESULTS

<u>SERA ID</u>	<u>1</u>	<u>2</u>	<u>MEAN</u>	<u># PASS</u>
15221	14.3	14.4	14.35	2/2
15242	41.3	37.5	39.4	2/2

COMMENTS: VISUAL CHECK OK.

APPROVED BY: *J. Ohta* DATE: 01/08/2009
 QA/QC