



**CERTIFICATE OF ANALYSIS**  
**(FINISHED KIT QUALITY CONTROL)**

PRODUCT NAME: **UREA NITROGEN (BUN) KINETIC**      MANUFACTURE DATE: 12/16/2009  
CATALOG #: B550-400  
LOT #: 22998      EXP. DATE: 8/2011

COMPONENTS:      1. BUN REAGENT      LOT#: 18420      EXP. DATE: 8/2011  
                         2. BUN STANDARD      LOT#: 19228      EXP. DATE: 10/2011

REAGENT BLANK VS. WATER: 1.525      AT 340 NM

**TOLERANCE LEVELS OF CONTROLS**

<u>SERA ID</u>	<u>EXP. DATE</u>	<u>TOLERANCE RANGE</u> (On TC-84)	<u>TARGET</u> (mg/dl)
15271	3/2012	11.2 – 16.8	14.0
15272	3/2012	36.4 – 54.6	45.5

**TEST VALUES/RESULTS**

<u>SERA ID</u>	<u>1</u>	<u>2</u>	<u>MEAN</u>	<u># PASS</u>
15271	12.3	14.6	13.45	2/2
15272	41.3	45.4	43.4	2/2

COMMENTS: VISUAL CHECK OK.

APPROVED BY:

  
Quality Control

DATE: 12/16/2009