

TECO DIAGNOSTICS



1268 N. Lakeview Ave., Anaheim, CA 92807, U.S.A.
Phone: (714) 463-1111 Toll Free (USA): (800) 222-9880
Fax: (714) 463-1169
www.tecodiag.com

CERTIFICATE OF ANALYSIS
(FINISHED KIT QUALITY CONTROL)

PRODUCT NAME: UREA NITROGEN (BUN) KINETIC MANUFACTURE DATE: 08/18/2010
CATALOG #: B550-240
LOT #: 27125 EXP. DATE: 08/2011

COMPONENTS: 1. BUN REAGENT LOT#: 18420 EXP. DATE: 08/2011
2. BUN STANDARD LOT#: 23729 EXP. DATE: 07/2012

REAGENT BLANK VS. WATER: AT 340 NM

TOLERANCE LEVELS OF CONTROLS

<u>SERA ID</u>	<u>EXP. DATE</u>	<u>TOLERANCE RANGE</u> (On TC-84)	<u>TARGET</u> (mg/dl)
15271	03/2012	11.2-16.8	14.0
15272	03/2012	36.4-54.6	45.5

TEST VALUES/RESULTS

<u>SERA ID</u>	<u>1</u>	<u>2</u>	<u>MEAN</u>	<u># PASS</u>
15271	15.1	13.6	14.35	2/2
15272	41.7	40.2	40.95	2/2

COMMENTS: VISUAL CHECK OK.

APPROVED BY: *[Signature]* DATE: 08/18/2010
Quality Control