



CERTIFICATE OF ANALYSIS
(FINISHED KIT QUALITY CONTROL)

PRODUCT NAME: UREA NITROGEN (BUN) KINETIC MANUFACTURE DATE: 02/09/2011
 CATALOG #: B550-400
 LOT #: 29338 EXP. DATE: 12/2012

COMPONENTS: 1. BUN REAGENT LOT#: 29874 EXP. DATE: 01/2014
 2. BUN STANDARD LOT#: 24042 EXP. DATE: 12/2012

REAGENT BLANK VS. WATER: 1.767 AT 340 NM

TOLERANCE LEVELS OF CONTROLS

<u>SERA ID</u>	<u>EXP. DATE</u>	<u>TOLERANCE RANGE</u> (On TC-84)	<u>TARGET</u> (mg/dl)
15281	09/2012	13.0-21.0	17.0
15272	03/2012	46.9-57.3	52.1

TEST VALUES/RESULTS

<u>SERA ID</u>	<u>1</u>	<u>2</u>	<u>MEAN</u>	<u># PASS</u>
15281	16	17	16.5	2/2
15272	47	47	47	2/2

COMMENTS: VISUAL CHECK OK.

APPROVED
 BY:


Quality Control

DATE: 02/09/2011