



CERTIFICATE OF ANALYSIS
(FINISHED KIT QUALITY CONTROL)

PRODUCT NAME: UREA NITROGEN (BUN) KINETIC MANUFACTURE DATE: 05/24/2011
 CATALOG #: B550-400
 LOT #: 32713 EXP. DATE: 10/2013

COMPONENTS: 1. BUN REAGENT LOT#: 29874 EXP. DATE: 01/2014
 2. BUN STANDARD LOT#: 31309 EXP. DATE: 10/2013

REAGENT BLANK VS. WATER: 1.770 AT 340 NM

TOLERANCE LEVELS OF CONTROLS

<u>SERA ID</u>	<u>EXP. DATE</u>	<u>TOLERANCE RANGE</u> (On TC-84)	<u>TARGET</u> (mg/dl)
15281	09/2012	13-21	17.3
15282	09/2012	39.9-59.9	49

TEST VALUES/RESULTS

<u>SERA ID</u>	<u>1</u>	<u>2</u>	<u>MEAN</u>	<u># PASS</u>
15281	15.3	17.6	16.45	2/2
15282	43.6	43.3	43.45	2/2

COMMENTS: VISUAL CHECK OK.

APPROVED
 BY:


Quality Control

DATE: 05/24/2011