



CERTIFICATE OF ANALYSIS
(FINISHED KIT QUALITY CONTROL)

PRODUCT NAME: UREA NITROGEN (BUN) KINETIC MANUFACTURE DATE: 07/14/2011
CATALOG #: B550-240
LOT #: 33628 EXP. DATE: 10/2013

COMPONENTS: 1. BUN REAGENT LOT#: 29874 EXP. DATE: 01/2014
2. BUN STANDARD LOT#: 31309 EXP. DATE: 10/2013

REAGENT BLANK VS. WATER: 1.565 AT 340 NM

TOLERANCE LEVELS OF CONTROLS

<u>SERA ID</u>	<u>EXP. DATE</u>	<u>TOLERANCE RANGE</u> (On TC-Matrix)	<u>TARGET</u> (mg/dl)
15281	09/2012	13-21	17
15282	09/2012	41-64	53

TEST VALUES/RESULTS

<u>SERA ID</u>	<u>1</u>	<u>2</u>	<u>MEAN</u>	<u># PASS</u>
15271	16	17	16.5	2/2
15282	48	48	48	2/2

COMMENTS: VISUAL CHECK OK.

APPROVED
BY:


Quality Control

DATE: 07/14/2011