



Certificate of Analysis

Product: Alkaline Phosphatase (AMP)
Catalogue Number: BXC0184C
Batch Number: 171477
Manufacture Date: 2017/06
Expiry Date: 2018/11

This product has been tested by Quality Control and has met all final QC specifications. This product is released for packing and is suitable for use until the stated expiry date.

The following results were obtained:

Control Description	Lot Number	Target (U/l)	Range (U/l)	Result (U/l)
Fortress Normal	HSN026	100	100-140	107.4
Fortress Elevated	HSE026	407	346-468	384

QC PASS/FAIL: **PASS**

Products containing Human Source Material have been tested at the donor level for HbsAg and antibodies to HIV1, HIV2 and HCV and found to be non-reactive.

Signed: Date: 2017/06

