



CERTIFICATE OF ANALYSIS

Cat No: BXC0321K
Lot no. 161665
Manufactured: 2016/07
Expiry: 2019/01
Product: Calibration Serum

PARAMETER	QUALITY ASSURANCE
Appearance Lyophilised	White lyophilised plug
Appearance after reconstitution	Yellow liquid
Appearance after 6 hours of manufacturing @ 25°C – open tube	No flocculation – maintains clarity
pH (25°C)	7.0
Bio-burden	<2 CFU/bottle

HOMOGENEITY across 20 vials	QUALITY ASSURANCE
AST	1.21 %
ALT	1.09 %
CK nac	1.01 %
Albumin	1.30%
Glucose	1.11%
Sodium	1.22%

This product was produced according to ISO 13485 guidelines.

QC PASS/FAIL: PASS

Signed: *Catherine Caskey*

Date: 2016/04

