

## Certificate of Analysis

Product: Direct LDL Cholesterol  
 Catalogue Number: BXC0431E  
 Batch Number: 161695  
 Manufactured: 2016/07  
 Expiry: 2017/10

This product has been tested by Quality Control and has met all final QC specifications. This product is released for packing and is suitable for use until the stated expiry date.

This product was produced according to ISO 13485 & ISO 9001 guidelines.

The following results were obtained:

Control Description	Lot Number	Target (mmol/l)	Range (mmol/l)	Result (mmol/l)
Lipid Control Level 1	LPCL01	1.09	0.90 – 1.30	1.11
Lipid Control Level 3	LPCE01	1.54	1.28 – 1.84	1.51

QC PASS/FAIL: **PASS**

Products containing Human Source Material have been tested at the donor level for HbsAg and antibodies to HIV1, HIV2 and HCV and found to be non-reactive.

Signed: *Catherine Caskey*

Date: 2016/07



ISO 13485:2003  
Cert. No. GB04 / 63491



ISO 9001:2000  
Cert. No. GB04 / 63490