



CERTIFICATE OF ANALYSIS
(FINISHED KIT QUALITY CONTROL)

PRODUCT NAME: **CALCIUM** MANUFACTURE DATE: 11/11/2010
 CATALOG #: C503-480
 LOT #: 28687 EXP. DATE: 11/2012

<u>COMPONENTS:</u>	<u>CATALOG #:</u>	<u>LOT#:</u>	<u>EXP. DATE:</u>
CALCIUM BUFFER REAGENT	C503-A	22775	12/2012
CALCIUM COLOR REAGENT	C503-B	21385	11/2012
CALCIUM STANDARD (10 mg/dl)	C503-C	21386	11/2012

BLANK: 0.210 AT 550 NM

TOLERANCE LEVELS OF SERUM CONTROLS

<u>SERA ID</u>	<u>EXP. DATE</u>	<u>TOLERANCE RANGE</u> (On TC-84)	<u>TARGET</u> (mg/dl)
15271	03/2012	7.2-8.8	8.0
15272	03/2012	10-12	11

TEST VALUES/RESULTS

<u>SERA ID</u>	<u>1</u>	<u>2</u>	<u>MEAN</u>	<u># PASS</u>
15271	8.5	7.9	8.2	2/2
15272	10.2	10.5	10.35	2/2

COMMENTS: VISUAL CHECK OK.

APPROVED BY: *V. Gascón* DATE: 11/11/201
 BY: Quality Control 0