



CERTIFICATE OF ANALYSIS
(FINISHED KIT QUALITY CONTROL)

PRODUCT NAME: **CALCIUM** MANUFACTURE DATE: 12/15/2010
 CATALOG #: C503-480
 LOT #: 29163 EXP. DATE: 11/2013

<u>COMPONENTS:</u>	<u>CATALOG #:</u>	<u>LOT#:</u>	<u>EXP. DATE:</u>
CALCIUM BUFFER REAGENT	C503-A	28816	11/2013
CALCIUM COLOR REAGENT	C503-B	28821	11/2013
CALCIUM STANDARD (10 mg/dl)	C503-C	28822	11/2013

BLANK: 0.748 AT 550 NM

TOLERANCE LEVELS OF SERUM CONTROLS

<u>SERA ID</u>	<u>EXP. DATE</u>	<u>TOLERANCE RANGE</u> (On TC-84)	<u>TARGET</u> (mg/dl)
15271	03/2012	7.2-8.8	8.0
15272	03/2012	10-12	11

TEST VALUES/RESULTS

<u>SERA ID</u>	<u>1</u>	<u>2</u>	<u>MEAN</u>	<u># PASS</u>
15271	7.7	7.5	7.6	2/2
15272	11.5	11.0	11.25	2/2

COMMENTS: VISUAL CHECK OK.

APPROVED BY: *V. Basay* DATE: 12/15/201
 BY: Quality Control 0