



**CERTIFICATE OF ANALYSIS**  
**(FINISHED KIT QUALITY CONTROL)**

PRODUCT NAME: **CHLORIDE** MANUFACTURE DATE: 04/12/2012  
CATALOG #: C501-480  
LOT #: 24797 EXP. DATE: 03/2013

COMPONENTS: CATALOG #: LOT#: EXP. DATE:  
CHLORIDE REAGENT C501-A 24800 03/2013  
CHLORIDE STANDARD (100 mEq/L) C501-B 24801 03/2013

BLANK: AT 505 NM

**TOLERANCE LEVELS OF SERUM CONTROLS**

<u>SERA ID</u>	<u>EXP. DATE</u>	<u>TOLERANCE RANGE</u>	<u>TARGET</u> (mEq/L)
15271	03/2012	77-110	99
15272	03/2012	78-102	90

**TEST VALUES/RESULTS**

<u>SERA ID</u>	<u>1</u>	<u>2</u>	<u>MEAN</u>	<u>#PASS</u>
15271	90.3	98.2	94.3	2/2
15272	87.0	84.8	85.9	2/2

COMMENTS: VISUAL CHECK OK.

APPROVED BY: *V. Passey* DATE: 04/12/2010  
Quality Control