



CERTIFICATE OF ANALYSIS
(FINISHED KIT QUALITY CONTROL)

PRODUCT NAME: **CHOLESTEROL LIQUID** MANUFACTURE DATE: 04/16/2010
 CATALOG #: C507-480
 LOT #: 25170 EXP. DATE: 09/2011

<u>COMPONENTS:</u>	<u>CATALOG #:</u>	<u>LOT#:</u>	<u>EXP. DATE:</u>
1. CHOLESTEROL LIQUID REAGENT	C507-A	25173	09/2011
2. CHOLESTEROL STANDARD (200 mg/dL)	C509-B	23346	12/2012

BLANK: 0.003 AT 505 NM

TOLERANCE LEVELS OF CONTROLS

<u>SERA ID</u>	<u>EXP. DATE</u>	<u>TOLERANCE RANGE</u>	<u>TARGET</u>
		(On TC-84)	(mg/dL)
15271	3/2012	221 – 290	245.5
15272	3/2012	90.7 – 110.9	100.8

TEST VALUES/RESULTS

<u>SERA ID</u>	<u>1</u>	<u>2</u>	<u>MEAN</u>	<u># PASS</u>
15271	260.5	243.4	252.0	2/2
15272	99.1	105.0	102.1	2/2

COMMENTS: VISUAL CHECK OK.

APPROVED BY:  DATE: 04/16/2010
 Quality Control