



**CERTIFICATE OF ANALYSIS**  
**(FINISHED KIT QUALITY CONTROL)**

PRODUCT NAME: **CHOLESTEROL LIQUID** MANUFACTURE DATE: 11/12/2010  
 CATALOG #: C507-480  
 LOT #: 28902 EXP. DATE: 04/2012

<u>COMPONENTS:</u>	<u>CATALOG #:</u>	<u>LOT#:</u>	<u>EXP. DATE:</u>
1. CHOLESTEROL LIQUID REAGENT	C507-A	28905	04/2012
2. CHOLESTEROL STANDARD (200 mg/dL)	C509-B	25174	04/2013

BLANK: 0.006 AT 505 NM

**TOLERANCE LEVELS OF CONTROLS**

<u>SERA ID</u>	<u>EXP. DATE</u>	<u>TOLERANCE RANGE</u> (On TC-84)	<u>TARGET</u> (mg/dL)
15271	03/2012	221(245.5)270	245.5
15272	03/2012	90.7(100.8)110.9	100.8

**TEST VALUES/RESULTS**

<u>SERA ID</u>	<u>1</u>	<u>2</u>	<u>MEAN</u>	<u># PASS</u>
15271	235.6	241	238.3	2/2
15272	106.8	94.2	100.5	2/2

COMMENTS: VISUAL CHECK OK.

APPROVED BY: *V. Basen* DATE: 11/12/2010  
 Quality Control