



CERTIFICATE OF ANALYSIS
(FINISHED KIT QUALITY CONTROL)

PRODUCT NAME: **CHOLESTEROL LIQUID** MANUFACTURE DATE: 05/23/2011
 CATALOG #: C507-480
 LOT #: 32744 EXP. DATE: 10/2012

<u>COMPONENTS:</u>	<u>CATALOG #:</u>	<u>LOT#:</u>	<u>EXP. DATE:</u>
1. CHOLESTEROL LIQUID REAGENT	C507-A	31308	10/2012
2. CHOLESTEROL STANDARD (200 mg/dL)	C509-B	30524	02/2014

BLANK: 0.003 AT 505 NM

TOLERANCE LEVELS OF CONTROLS

<u>SERA ID</u>	<u>EXP. DATE</u>	<u>TOLERANCE RANGE</u> (On TC-Matrix)	<u>TARGET</u> (mg/dL)
15281	09/2012	201-314	258
15282	09/2012	85-130	108

TEST VALUES/RESULTS

<u>SERA ID</u>	<u>1</u>	<u>2</u>	<u>MEAN</u>	<u># PASS</u>
15281	246	245	245.5	2/2
15282	103	102	102.5	2/2

COMMENTS: VISUAL CHECK OK.

APPROVED BY:  DATE: 05/23/2011
 Quality Control