



CERTIFICATE OF ANALYSIS
(FINISHED KIT QUALITY CONTROL)

PRODUCT NAME: **CHOLESTEROL LIQUID** MANUFACTURE DATE: 11/10/2011
 CATALOG #: C507-480
 LOT #: 35784 EXP. DATE: 04/2013

| <u>COMPONENTS:</u> | <u>CATALOG #:</u> | <u>LOT#:</u> | <u>EXP. DATE:</u> |
|-------------------------------------|-------------------|--------------|-------------------|
| 1. CHOLESTEROL LIQUID REAGENT | C507-A | 35783 | 04/2013 |
| 2. CHOLESTEROL STANDARD (200 mg/dL) | C509-B | 34814 | 09/2014 |

BLANK: 0.017 AT 505 NM


TOLERANCE LEVELS OF CONTROLS

| <u>SERA ID</u> | <u>EXP. DATE</u> | <u>TOLERANCE RANGE</u> | <u>TARGET</u> (mg/dL) |
|----------------|------------------|------------------------|--------------------------|
| 15281 | 09/2012 | 201 - 314 | 258 |
| 15282 | 09/2012 | 85 - 130 | 108 |

TEST VALUES/RESULTS

| <u>SERA ID</u> | <u>1</u> | <u>2</u> | <u>MEAN</u> | <u># PASS</u> |
|----------------|----------|----------|-------------|---------------|
| 15281 | 257 | 254 | 256 | 2/2 |
| 15282 | 106 | 105 | 106 | 2/2 |

COMMENTS: VISUAL CHECK OK.

APPROVED BY:  DATE: 11/10/2011
 Nour Qasqas
 Quality Control Lab. Director