

TECO DIAGNOSTICS



1268 N. Lakeview Ave., Anaheim, CA 92807, U.S.A.
Phone: (714) 463-1111 Toll Free (USA): (800) 222-9880
Fax: (714) 463-1169
www.tecodiag.com

CERTIFICATE OF ANALYSIS
(FINISHED KIT QUALITY CONTROL)

PRODUCT NAME: **CHOLESTEROL - 10** MANUFACTURE DATE: 08/26/2010
CATALOG #: C509-150
LOT #: 27298 EXP. DATE: 04/2013

COMPONENTS: 1.CHOLESTEROL REAGENT LOT #: 25353 EXP. DATE: 05/2013
2.CHOLESTEROL STANDARD LOT #: 25174 EXP. DATE: 04/2013

BLANK: 0.015 AT 505 NM

STANDARD: 200 mg/dL


TOLERANCE LEVELS OF SERUM CONTROL

| <u>SERA ID</u> | <u>EXP. DATE</u> | <u>TOLERANCE RANGE</u> (On TC-MATRIX) | <u>TARGET</u> (mg/dL) |
|----------------|------------------|--|--------------------------|
| 15271 | 02/2009 | 217.8 – 266.2 | 242 |
| 15272 | 02/2009 | 90.5 – 110.7 | 100.6 |

TEST VALUES/RESULTS

| <u>SERA ID</u> | <u>1</u> | <u>2</u> | <u>MEAN</u> | <u># PASS</u> |
|----------------|----------|----------|-------------|---------------|
| 15271 | 242 | 240 | 241 | 2/2 |
| 15272 | 103 | 104 | 103.5 | 2/2 |

COMMENTS: VISUAL CHECK OK.

| | | | |
|------------------|--|-----------|-------------------|
| APPROVE D BY: |  Quality Control | DATE : | <u>08/26/2010</u> |
|------------------|--|-----------|-------------------|