

TECO DIAGNOSTICS



1268 N. Lakeview Ave., Anaheim, CA 92807, U.S.A.  
Phone: (714) 463-1111 Toll Free (USA): (800) 222-9880  
Fax: (714) 463-1169  
[www.tecodiag.com](http://www.tecodiag.com)

**CERTIFICATE OF ANALYSIS**  
**(FINISHED KIT QUALITY CONTROL)**

PRODUCT NAME: **CHOLINESTERASE** MANUFACTURE DATE: 05/28/2010  
CATALOG #: C511-60  
LOT #: 25834 EXP. DATE: 08/2012

COMPONENT: CHOLINESTERASE REAGENT LOT#: 25481 EXP. DATE: 08/2012

REAGENT BLANK VS. WATER: 0.236 AT 405 NM

**TOLERANCE LEVELS OF SERUM CONTROLS**

<u>SERA ID</u>	<u>EXP. DATE</u>	<u>TOLERANCE RANGE</u>	<u>TARGET (IU/L)</u>
15271	03/2012	1448-1767	1607
15272	03/2012	1751-2140	1946

**TEST VALUES/RESULTS**

<u>SERA ID</u>	<u>1</u>	<u>2</u>	<u>MEAN</u>	<u># PASS</u>
15271	1507	1551	1529	2/2
15272	1798	1904	1851	2/2

COMMENTS: VISUAL CHECK OK.

APPROVED BY:  DATE: 05/28/2010  
Quality Control