



CERTIFICATE OF ANALYSIS
(FINISHED KIT QUALITY CONTROL)

PRODUCT NAME: **CHOLINESTERASE** MANUFACTURE DATE: 08/01/2011
CATALOG #: C511-60
LOT #: 33891 EXP. DATE: 04/2013

COMPONENT: CHOLINESTERASE REAGENT LOT#: 29689 EXP. DATE: 04/2013

REAGENT BLANK VS. WATER: 0.216 AT 405 NM

TOLERANCE LEVELS OF SERUM CONTROLS

<u>SERA ID</u>	<u>EXP. DATE</u>	<u>TOLERANCE RANGE</u>	<u>TARGET (IU/L)</u>
15281	09/2012	2988-4432	3735
15282	09/2012	4713-7069	5891

TEST VALUES/RESULTS

<u>SERA ID</u>	<u>1</u>	<u>2</u>	<u>MEAN</u>	<u># PASS</u>
15281	4037	4027	4032	2/2
15282	6476	6233	6354	2/2

COMMENTS: VISUAL CHECK OK.

APPROVED BY: 
Quality Control DATE: 08/03/2011