



CERTIFICATE OF ANALYSIS
(FINISHED KIT QUALITY CONTROL)

PRODUCT NAME: **CK – NAC (KINETIC)** MANUFACTURE DATE: 11/09/2010
 CATALOG #: C512-60
 LOT #: 28691 EXP. DATE: 08/2012

COMPONENT: CK-NAC REAGENT (CAT#: C512-A) LOT#: 23363 EXP. DATE: 08/2012

REAGENT BLANK VS. WATER: 0.349 AT 340 NM

TOLERANCE LEVELS OF SERUM CONTROLS

<u>SERA ID</u>	<u>EXP. DATE</u>	<u>TOLERANCE RANGE</u>	<u>TARGET</u> (IU/L)
15271	03/2012	105.5 – 158.2	131.9
15272	03/2012	371.9 – 557.9	464.9

TEST VALUES/RESULTS

<u>SERA ID</u>	<u>1</u>	<u>2</u>	<u>MEAN</u>	<u># PASS</u>
15271	146.2	142.5	144.35	2/2
15272	490.5	482.2	486.35	2/2

COMMENTS: VISUAL CHECK OK.

APPROVED
BY:


Quality Control

DATE: 11/09/2010