



CERTIFICATE OF ANALYSIS
(FINISHED KIT QUALITY CONTROL)

PRODUCT NAME: **CK-MB REAGENT SET** MANUFACTURE DATE: 05/27/2010
 CATALOG #: C614-60
 LOT #: 25956 EXP. DATE: 06/2012

COMPONENTS: 1. CK-MB REAGENT LOT#: 23363 EXP. DATE: 08/2012
 2. CK-MB DILUENT LOT#: 25959 EXP. DATE: 06/2012

TOLERANCE LEVELS OF SERUM CONTROLS

<u>CONTROLS</u>	<u>LOT#</u>	<u>EXP. DATE</u>	<u>TOLERANCE RANGE</u> (TC MATRIX)	<u>TARGET</u> (IU/L)
LEVEL I	I-68-31	12/2010	<25	11.5
LEVEL II	I-72-5	08/2011	40-60	50
LEVEL III	I-72-6	03/2012	73-109	91

TEST VALUES/RESULTS

<u>SERA ID</u>	<u>1</u>	<u>2</u>	<u>MEAN</u>	<u># PASS</u>
LEVEL I	89	92	91	2/2
LEVEL II	48	48	48	2/2
LEVEL III	4	4	4	2/2

COMMENTS: VISUAL CHECK OK.

APPROVED
BY:


Quality Control

DATE: 05/27/2010