

TECO DIAGNOSTICS

1268 N. Lakeview Ave., Anaheim, CA 92807, U.S.A.

Phone: (714) 693-7788 Toll Free: 1-800-222-9880

Fax: (714) 693-3838

www.tecodiag.com Email: tecodiag@tecodiag.com**CERTIFICATE OF ANALYSIS**
(FINISHED KIT QUALITY CONTROL)PRODUCT NAME: **CK-MB**
CATALOG #: C614-30
LOT #: 34639MANUFACTURE DATE: 08/31/2011
EXP. DATE: 06/2013

COMPONENTS:	CK-MB REAGENT	LOT#:	30901	EXP. DATE:	09/2013
	CK-MB DILUENT	LOT#:	32751	EXP. DATE:	06/2013

TOLERANCE LEVELS OF SERUM CONTROLS

<u>CONTROLS</u>	<u>LOT#</u>	<u>EXP. DATE</u>	<u>TOLERANCE RANGE</u> (TC MATRIX)	<u>TARGET</u>
LEVEL I	I-74-25	12/2012	<25	11.5
LEVEL II	I-74-26	12/2012	40 – 60	50
LEVEL III	I-72-6	03/2012	73 – 109	91

TEST VALUES/RESULTS

<u>SERA ID</u>	<u>1</u>	<u>2</u>	<u>MEAN</u>	<u># PASS</u>
LEVEL I	16	18	17	2/2
LEVEL II	43	41	42	2/2
LEVEL III	106	106	106	2/2

COMMENTS: VISUAL CHECK OK

APPROVED BY:


Quality Control
DATE: 08/31/2011