



CERTIFICATE OF ANALYSIS
(FINISHED KIT QUALITY CONTROL)

PRODUCT NAME: **CREATININE (KINETIC)** MANUFACTURE DATE: 03/23/2010
 CATALOG #: C514-480
 LOT #: 24820 EXP. DATE: 06/2012

<u>COMPONENTS:</u>	<u>CATALOG #:</u>	<u>LOT #:</u>	<u>EXP. DATE:</u>
1. PICRIC ACID REAGENT	C514-A	24707	03/2013
2. SODIUM HYDROXIDE (NaOH) REAGENT	C514-B	24649	03/2013
3. STANDARD (5.0 mg/dL)	C514-C	23443	06/2012

BLANK: 0.119 AT 505 NM

TOLERANCE LEVELS OF SERUM CONTROL

<u>SERA ID</u>	<u>EXP. DATE</u>	<u>TOLERANCE RANGE</u> (On TC-Matrix)	<u>TARGET</u> (mg/dL)
15271	03/2012	0.9-1.5	1.2
15272	03/2012	4.9-6.5	5.6

TEST VALUES/DATA

<u>SERA ID</u>	<u>1</u>	<u>2</u>	<u>MEAN</u>	<u># PASS</u>
15271	1.15	1.22	1.19	2/2
15272	5.19	5.89	5.54	2/2

COMMENTS: VISUAL CHECK OK.

APPROVED BY:  DATE: 03/23/2010
 Quality Control