
CERTIFICATE OF ANALYSIS
(FINISHED KIT QUALITY CONTROL)

PRODUCT NAME: **CREATININE (KINETIC)** MANUFACTURE DATE: 09/07/2011
CATALOG #: C514-480
LOT #: 34526 EXP. DATE: 11/2013

<u>COMPONENTS:</u>	<u>CATALOG #:</u>	<u>LOT #:</u>	<u>EXP. DATE:</u>
1. PICRIC ACID REAGENT	C514-A	33404	07/2014
2. SODIUM HYDROXIDE (NaOH) REAGENT	C514-B	33405	07/2014
3. STANDARD (5.0 mg/dL)	C514-C	33016	11/2013

BLANK: 0.154 AT 505 NM

TOLERANCE LEVELS OF SERUM CONTROL

<u>SERA ID</u>	<u>EXP. DATE</u>	<u>TOLERANCE RANGE</u>	<u>TARGET</u> (mg/dL)
15281	09/2012	1.0 – 1.5	1.3
15282	09/2012	4.5 – 6.7	5.6

TEST VALUES/DATA

<u>SERA ID</u>	<u>1</u>	<u>2</u>	<u>MEAN</u>	<u># PASS</u>
15281	1.2	1.3	1.3	2/2
15282	5.8	5.8	5.8	2/2

COMMENTS: VISUAL CHECK OK.

APPROVED BY:  DATE: 09/07/2011
Quality Control