

TECO DIAGNOSTICS



1268 N. Lakeview Ave., Anaheim, CA 92807, U.S.A.
Phone: (714) 463-1111 Toll Free (USA): (800) 222-9880
Fax: (714) 463-1169
www.tecodiag.com

CERTIFICATE OF ANALYSIS
(FINISHED KIT QUALITY CONTROL)

PRODUCT NAME: GAMMA-GT (SOLUBLE) MANUFACTURE DATE: 05/21/2010
CATALOG #: G571-150
LOT#: 25594 EXP. DATE : 07/2012

COMPONENTS: GGT SOLUBLE LOT# 23999 EXP. DATE: 07/2012

WORKING REAGENT VS. DI WATER: 0.380 AT 405 NM

TOLERANCE LEVELS OF SERUM CONTROL

<u>SERA ID</u>	<u>EXP. DATE</u>	<u>TOLERANCE RANGE</u>	<u>TARGET</u> (IU/L)
15271	03/2012	23.7 – 35.5	29.6
15272	03/2012	96.1 – 172	134.1

TEST DATA

<u>SERA ID</u>	<u>1</u>	<u>2</u>	<u>MEAN</u>	<u># PASS</u>
15271	28.8	28.3	28.6	2/2
15272	107.2	109.5	108.4	2/2

COMMENTS: VISUAL CHECK OK.

APPROVED BY:


Quality Control

DATE: 05/21/2010