

TECO DIAGNOSTICS



1268 N. Lakeview Ave., Anaheim, CA 92807, U.S.A.  
Phone: (714) 463-1111 Toll Free (USA): (800) 222-9880  
Fax: (714) 463-1169  
[www.tecodiag.com](http://www.tecodiag.com)

**CERTIFICATE OF ANALYSIS**  
**(FINISHED KIT QUALITY CONTROL)**

PRODUCT NAME: GAMMA-GT (SOLUBLE) MANUFACTURE DATE: 11/17/2010  
CATALOG #: G571-204  
LOT#: 29075 EXP. DATE : 03/2013

COMPONENTS: GGT SOLUBLE LOT# 27061 EXP. DATE: 03/2013

WORKING REAGENT VS. DI WATER: 0.392 AT 405 NM

**TOLERANCE LEVELS OF SERUM CONTROL**

<u>SERA ID</u>	<u>EXP. DATE</u>	<u>TOLERANCE RANGE</u>	<u>TARGET</u> (IU/L)
15271	03/2012	23.7(29.6)35.5	29.6
15272	03/2012	96.1(134.1)172.1	134.1

**TEST DATA**

<u>SERA ID</u>	<u>1</u>	<u>2</u>	<u>MEAN</u>	<u># PASS</u>
15271	34.4	33.9	34.15	2/2
15272	111.1	117.4	114.25	2/2

COMMENTS: VISUAL CHECK OK.

APPROVED  
BY:

  
Quality Control

DATE: 11/17/2010