

TECO DIAGNOSTICS



1268 N. Lakeview Ave., Anaheim, CA 92807, U.S.A.
Phone: (714) 463-1111 Toll Free (USA): (800) 222-9880
Fax: (714) 463-1169
www.tecodiag.com

CERTIFICATE OF ANALYSIS
(FINISHED KIT QUALITY CONTROL)

PRODUCT NAME: GAMMA-GT (SOLUBLE) MANUFACTURE DATE: 07/22/2014
CATALOG #: G571-150
LOT#: 56271 EXP. DATE: 11/2016

COMPONENTS: GGT SOLUBLE LOT# 54708 EXP. DATE: 11/2016

TOLERANCE LEVELS OF SERUM CONTROL

<u>SERA ID</u>	<u>EXP. DATE</u>	<u>TOLERANCE RANGE</u> (on TC-Matrix)	<u>TARGET</u> (IU/L)
27611	05/2015	26.9 – 40.3	33.6
27602	03/2015	109.2 – 163.8	136.5

TEST DATA

<u>SERA ID</u>	<u>MEAN</u>	<u>PASS/FAIL</u>
27611	31.4	PASS
27602	121.0	PASS

COMMENTS: VISUAL CHECK OK.

APPROVED BY: Susan Tang
Susan Tang
Quality Control Manager

DATE: 07/23/2014