



CERTIFICATE OF ANALYSIS
(FINISHED KIT QUALITY CONTROL)

PRODUCT NAME: **GLUCOSE OXIDASE** MANUFACTURE DATE: 05/27/2010
 CATALOG #: G519-1L
 LOT #: 25953 EXP. DATE: 05/2012

<u>COMPONENTS</u>	<u>CATALOG #</u>	<u>LOT #</u>	<u>EXP. DATE</u>
GLUCOSE OXIDASE REAGENT	G519-A	19420	05/2012
GLUCOSE STANDARD (100 mg/dL)	G518-B	21637	11/2012

BLANK: -0.010 AT 505 NM

TOLERANCE LEVELS OF CONTROLS

<u>SERA ID</u>	<u>EXP. DATE</u>	<u>TOLERANCE RANGE</u>	<u>TARGET</u> (mg/dL)
15271	03/2012	71 - 107	89
15272	03/2012	237 - 355	296

TEST VALUES/RESULTS

<u>SERA ID</u>	<u>1</u>	<u>2</u>	<u>MEAN</u>	<u># PASS</u>
15271	9036	106.9	98.8	2/2
15272	317.8	274.4	296.1	2/2

COMMENTS: VISUAL CHECK OK.

APPROVE D BY:	 Quality Control	DATE :	<u>05/27/2010</u>
------------------	--	-----------	-------------------