



**CERTIFICATE OF ANALYSIS**  
**(FINISHED KIT QUALITY CONTROL)**

PRODUCT NAME: **GLUCOSE OXIDASE**      MANUFACTURE DATE: 11/10/2011  
 CATALOG #: G519-1L  
 LOT #: 35787      EXP. DATE: 06/2014

<u>COMPONENTS</u>	<u>CATALOG #</u>	<u>LOT #</u>	<u>EXP. DATE</u>
GLUCOSE OXIDASE REAGENT	G519-A	32179	06/2014
GLUCOSE STANDARD (100 mg/dL)	G518-B	34317	08/2014

BLANK: 0.012 AT 505 NM


**TOLERANCE LEVELS OF CONTROLS**

<u>SERA ID</u>	<u>EXP. DATE</u>	<u>TOLERANCE RANGE</u>	<u>TARGET</u> (mg/dL)
15281	09/2012	71 - 112	92
15282	09/2012	229 - 377	303

**TEST VALUES/RESULTS**

<u>SERA ID</u>	<u>1</u>	<u>2</u>	<u>MEAN</u>	<u># PASS</u>
15281	87	87	87	2/2
15282	277	279	278	2/2

COMMENTS: VISUAL CHECK OK.

APPROVED BY:  DATE: 11/10/2011  
 Nour Qasqas  
 Quality Control Lab. Director